

Experience **Simple Connected Health**.

To address the needs of home health patients, data flows instantly into existing workflows without the need for pairing, Wi-Fi, or apps resulting in a frictionless, integrated workflow for both patients and care teams.

Stay connected using the **Coag-Sense® PT/INR Meter combined with the Stel Home Vitals Hub**. To set-up your connection, follow these easy step-by-step instructions:

1. Plug in your Stel Vitals Hub and wait for the blinking red light to become solid.
2. Power On the fully charged PT/INR meter by pressing and holding the power button on the side.
3. Using the touchscreen, press the 'Settings' button (⚙️) in the bottom right corner.
4. Within the settings menu, press 'Communication settings.'
5. In the Communication settings menu, press 'Bluetooth.' Two options will appear, 'NO RESET' and 'RESET.' To pair, press the 'RESET' button and wait for the passcode screen to appear.
6. Using the number pad, enter the passcode '573171' then press the right arrow '>' to submit.
7. After a few moments, the INR will tell you know it has connected to the Hub. If pairing fails, return to the home page using the square button in the center at the bottom of the screen and repeat steps 2 - 6. If pairing succeeds, press the square button in the center at the bottom of the screen to return to the HOME page and take a measurement!
8. If you are away from your Stel Vitals Hub:
 - a. Take your INR measurements at the frequency recommended by your care provider.
 - b. When you are back in the range of the Hub, simply take your next INR measurement.
 - c. When you see the BLUE light blink, all the readings taken while you were away from your Stel Vitals Hub will be transmitted.

If you have any questions about the Stel Hub or if you don't see a blinking blue light after taking a measurement, **please contact your Clinical Care Team**.