



AEGISPOC 4.17.1 CUSTOMER RELEASE NOTES

DEVICES, UPDATES, AND CORRECTIONS

The release of AegisPOC™ 4.17.1 includes a device update.
For a summary of all updates, please reference the list
below.



DEVICES

- Added new device driver:
 - Vendor - CoaguSense, Inc.
 - Device - Coag-Sense®

DEVICE DETAILS

- The Coag-Sense Prothrombin Time (PT)/INR Monitoring System is an in vitro diagnostic device that provides quantitative prothrombin time (PT) results, expressed in seconds and International Normalized Ratio (INR) units.
- It uses fresh capillary whole blood. It is intended for use by health care professionals at the point of care to monitor patients who are on warfarin-type (coumarin) anticoagulation therapy. The device is not intended to be used for screening purposes

TO LEARN MORE, VISIT [GLOBALPOINTOFCARE.ABBOTT](https://globalpointofcare.abbott)

Products not available in all countries.

© 2024 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.
Any photos displayed are for illustrative purposes only. AI-DSREL0995 Rev. 1.0 09/2024

